

Service manual

Refrigerator/freezer

BL-Series PRO-ACTIVE / ULTIMATE



Read the instructions prior to performing any task!

Philipp Kirsch GmbH
Im Lossenfeld 14
77731 Willstätt-Sand
GERMANY
Telephone: +49 781 9227-0
Fax: +49 781 9227-200
Email: info@kirsch-medical.de
Internet: www.kirsch-medical.de

KD_Typ_XXX_ab Fabriknummer_XXX XX XXXXX, 1, en_GB

Target audience

This manual is aimed at the following persons:

- Personnel:
- Customer service
 - Certified refrigeration technician

Models

This manual is valid for the following model:

Area	Model	Serial number from
Blood cooling	BL 720 PRO-ACTIVE	700 31 25002

The unit has the following specifications:

- independent rapid defrosting using hot gas
- Forced-air cooling
- Melt water evaporation
- Cold room lighting (optional)
- Delayed inner temperature display
- Lockable door
- Key switch for control
- Door handle
- Self-closing door
- Optical and acoustic warning system
- Power failure warning

Contact

Germany

PHILIPP KIRSCH GMBH

Im Lossenfeld 14

77731 Willstätt-Sand, Germany

Telephone: +49 (0) 781 9227 0

Fax: +49 (0) 781 9227 200

Email: info@kirsch-medical.de

Mon - Thurs: 08.00 to 12.15, 13.15 to 16.30 Fri: 08.00 to 12.00, 13.00 to 16.00

Netherlands

VDW CoolSystems

De Aaldor 8c | NL-4191 PC Geldermalsen, Nederland

Telephone: +31 345 582 546

Fax: +31 345 580 390

E-mail: service@vdwcoolsystems.nl

Order acceptance: 24/7 (24 hours, seven days a week)

Illustrations

The illustrations in this manual are examples. They are intended for a basic understanding and may vary depending on the model.

Table of contents

1	Safety	7
1.1	Symbols in this manual	7
1.2	Intended use	8
1.3	Foreseeable misuse	8
1.4	Residual risks	8
1.5	Safety labelling	11
1.6	Personnel qualification	12
1.7	Personal protective equipment	13
2	Detailed function description	15
2.1	Details on control	15
2.2	Temperature sensor	17
2.3	Error message and buzzer	17
2.4	Safety regulation	18
2.5	Key switch	19
2.6	Defrosting	19
2.7	Data export on the device with PC-KIT-STICK	20
3	Function check	23
3.1	Checking the alarm trigger	23
3.2	Checking sensors	23
3.3	Initialising the control and monitoring unit	24
3.4	Comparison of actual values	24
4	Maintenance and Customization	27
4.1	Status displays and error messages	27
4.1.1	Status displays	27
4.1.2	Error messages	28
4.1.3	Setting the time and date	33
4.2	Switching the door hinge	34
4.3	Switching the door handle	35
4.4	Replacing the display unit	35
4.5	Replacing sensors	35
4.6	Replacing the control and monitoring unit	36
4.7	Replacing the interior blower	36
4.8	Replacing the door lock	36
4.9	Replacing the door seal	37
4.10	Replacing the door contact	37
4.11	Replacing the battery	37
4.12	Replacing the cold room lighting (optional)	37
4.13	Replacing the cooling machine	38
5	Customizing parameters	41
6	Appendix	45
6.1	Technical data	46
6.2	wiring diagram	49

6.3	Exploded view.....	51
6.4	Spare parts list.....	52
7	Index.....	57

1 Safety

This section provides an overview of all important safety aspects for optimal protection of patients and staff, and for safe and trouble-free operation of the unit.

Non-compliance with the instructions and safety warnings in this instruction manual can cause considerable risks.

1.1 Symbols in this manual

Safety instructions

Safety instructions are marked with symbols in this manual. Safety instructions are initiated by signal words that express the degree of risk.

In order to avoid accidents, injury or damage and ensure maximum patient safety, always comply with safety instructions and act with care.



DANGER!

This combination of symbol and signal word indicates an imminently hazardous situation that will lead to death or serious injury unless avoided.



WARNING!

This combination of symbol and signal word indicates a potentially hazardous situation that can lead to death or serious injury unless avoided.



CAUTION!

This combination of symbol and signal word indicates a potentially hazardous situation that can lead to minor or slight injury unless avoided.



NOTICE!





This combination of symbol and signal word indicates a potentially hazardous situation that can lead to property damage or environmental damage unless avoided.

Hints and recommendations



This symbol highlights useful hints and recommendations as well as information for efficient and trouble-free use of the unit.

Other markings

Mark	Explanation
	Step-by-step instructions
	Results of actions
	References to sections in this manual
	Lists without a specified order

1.2 Intended use

The blood storage refrigerator is used for the commercial storage of blood reserves. The blood bank refrigerator is a medical product.

1.3 Foreseeable misuse

The unit is not designed for domestic use. The unit is used for commercial storage of chilled goods in line with its intended purpose.

Do not use the unit to cool warm goods. Do not store chilled goods in the unit if their cooling chain was interrupted during delivery or stock transfer.

Do not store food or drink in the unit.

Do not store chilled goods that exceed the carrying capacity of the wire shelves and drawers.

1.4 Residual risks

Hot surface



WARNING!

Danger due to hot surface!

The marked areas of the unit can cause severe skin injuries if touched.

- Do not touch areas of the unit marked in this way.
- These areas are very hot and can still cause burns several hours after the unit has been switched off.

Blocked interior ventilation



WARNING!

Risk of chilled goods damage due to a blocked inner ventilation !

A blocked inner ventilation can cause issues to the air circulation inside the unit so that temperature homogeneity won't be ensured any longer. This can lead to a massive temperature deviation inside the unit and may result in damage to the chilled goods.

- Make sure that the ventilation in the upper area is not blocked.
- Do not cover the ventilation grids with chilled goods.
- Make sure that the air outlet in the lower area is not blocked.

Escaping refrigerant



WARNING!

Danger due to escaping refrigerant!

The refrigerant used in the unit is explosive and pressurised, and can cause serious injury if it comes into contact with the eyes or skin.

- During transport and set-up, do not bend or pierce the tubing and the evaporator.
- Do not damage the surface coat (scratching it off, for example).
- Wear safety glasses and protective gloves when handling the refrigerant circuit.

Tilting the unit



WARNING!

Danger of crushing injuries from falling unit!

If the unit is tilted, it will fall over in an uncontrolled fashion. When the unit falls over, there is a danger of crushing to the hands and feet.

- Wear safety shoes and safety gloves when transporting the unit.
- Transport the unit in an upright position.

Missing power supply



CAUTION!

Danger of goods and personal damage due to missing power supply!

In the event of a power failure, the battery supplies power only to the warning system.

The cooling unit is switched off, which can damage or destroy the chilled goods and lead to a serious injuries as a result

- After a power failure, make sure that the cooling unit is supplied with power again.
- If necessary, move the chilled goods.

Contamination of chilled goods



CAUTION!

Danger to chilled goods due to contamination!

A dirty unit can contaminate the stored chilled goods. A contamination can damage or destroy the chilled goods.

- Do not store contaminated chilled goods in the unit.
- After any contamination, clean, disinfect and sterilise the unit and the chilled goods.
- Clean, disinfect and sterilise the unit in line with this instruction manual during operation.
- Comply with the applicable requirements for personal hygiene, disinfection and sterilisation.

Frostbite from refrigerant



WARNING!

Risk of frostbite from refrigerant!


At normal atmospheric pressure and ambient temperatures, liquid refrigerant vaporises so quickly that contact with skin or eyes could cause frostbite (risk of blinding).





- Wear safety glasses and protective gloves when handling the refrigerant circuit.
- Never carry out work on the cooling circuit independently. Never force the cooling circuit open.
- Please pay attention to the type plates in the interior and on the cooling machine.

1.5 Safety labelling

Rating plate

Modell 1 Type	BL 100 PRO-ACTIVE	Kühl-Inhalt Gross volume	95 Liter	6
2 SN	100 32 35000	UDI	42606885800101003235000	7
3 Kältemittel Refrigerant	R600a 0,03 kg	Kühl-Aggregat Cooling system	VKD 2713 K	8
4 Wechselstrom Voltage	220-240V~50 Hz	Stromaufnahme Power input	0,31 A	9
5 Nutzraumtemperaturbereich Usable- space temperature	+4°C	Umgebungstemperaturbereich Ambient temperature	10°bis 38°C	10











Philipp Kirsch GmbH
 Im Lössenfeld 14
 DE-77731, Willstätt-Sand
 Tel.: +49 781 9227-0
 Web: www.kirsch-medical.com

Medical device class IIa
 Regulation (EU) 2017/745

Fig. 1: Example of a rating plate (BL 100 PRO-ACTIVE)

The rating plate contains the following information (information corresponds to the red number):

- 1 - Model
- 2 - Serial number
- 3 - Refrigerant
- 4 - Alternating current
- 5 - Usable room temperature range
- 6 - Cooling system capacity
- 7 - UDI
- 8 - Refrigeration unit
- 9 - Power consumption
- 10 - Ambient temperature range

Symbol	Meaning
	CE labelling
	Read the instructions for use
	Manufacturer
	Date of manufacture
	Is a medical device
	Serial number
	Attention
	Unique product identification



Old devices must be disposed of separately from household waste.

Batteries or rechargeable batteries that are not permanently installed must be removed beforehand and disposed of separately.

1.6 Personnel qualification

Insufficiently qualified personnel



WARNING!

Danger of damage and injury due to commissioning by unqualified staff!

If commissioning is performed incorrectly by unqualified staff, serious damage to the chilled goods can result, which in turn can seriously harm patients.

- Have all tasks performed only by staff qualified for those tasks.
- Keep unauthorised persons away from the working area.

Personnel qualifications

Medical devices may only be set up, operated, used and maintained by persons who have the necessary training or knowledge and experience.

These instructions specify the personnel qualifications required for the various areas of activity listed below:

Certified refrigeration technician

Due to their specialist qualification, knowledge and experience, and knowledge of the applicable standards and regulations, certified refrigeration technicians are able to perform work on refrigeration systems and the refrigerant in question, and identify and avoid potential dangers of their own accord.

The acquired certification also includes the expertise required to prevent emissions and safely handle refrigeration equipment of the relevant type and size. Refrigeration technicians have been trained and are certified to work at the specific location of use and are familiar with the relevant standards and regulations.

Customer service

Certain work can only be performed by the manufacturer's customer service or by a specialist company contracted by the manufacturer. No other staff is permitted to carry out this work. Contact the manufacturer's customer service to carry out the necessary work.

Kirsch service contact:

PHILIPP KIRSCH GMBH

Im Lössenfeld 14

77731 Willstätt-Sand, Germany

Telephone: +49 (0) 781 9227 777

Fax: +49 (0) 781 9227 200

Email: info@kirsch-medical.de

Mon – Thurs: 8.00am to 12.15pm, 1.15pm to 4.30pm, Fri: 8.00am to 12.00pm, 1.00pm to 4.00pm

General personnel qualification requirements

Only persons who can be expected to carry out their work in a reliable manner are authorised as personnel. Persons whose ability to react is affected by, for example, drugs, alcohol or medication are not authorised.

When selecting personnel, please take into account the age and job-related regulations in force at the place of work.

1.7 Personal protective equipment

Personal protective equipment is used to protect personnel against hazards that could affect their safety or health at work.

Always put on the personal protective equipment required in the various chapters of these instructions before starting the work in question.

Follow the instructions on personal protective equipment in the work area.

Personnel must wear personal protective equipment when performing the various tasks on and with the device. This equipment is referred to separately in the individual chapters in these instructions. This personal protective equipment is explained below:



NOTICE!

Wear suitable protective gloves when using the unit. These gloves must comply with either PPE Ordinance 2016/424 or the MDD Directive 93/42/EEC (superseeded on 27 May 2021 by MDR Regulation 2017/745).



Chemical-resistant protective gloves

Chemical-resistant protective gloves protect the hands from aggressive chemicals.



Protective gloves

Protective gloves protect the hands from friction, abrasions, piercing or deeper injuries as well as from contact with hot surfaces.



Protective goggles

Protective goggles serve to protect the eyes from flying debris and splashes of liquid.

2 Detailed function description

2.1 Details on control




Fig. 2: Display and control unit

- A Key switch
- B Display and control unit (variable key assignment)
- C USB port










Please find the functions description of the buttons and signs listed below:

Tab. 1: Buttons

Pos. no.	Button	Description	Function
1		[Defrosting]	Activate additional defrosting.
		[SuperFrost]	Activate SuperFrost function.
		[Humidity]	Adjust humidity.
2		[Max.]	Showing maximum value of the temperature memory.
3		[Reset]	Reset temperature memory. Switch off buzzer.

Detailed function description




Details on control


Pos. no.	Button	Description	Function
4		<i>[Temperature warning max.]</i>	Set point upper temperature limit
5		<i>[Start]</i>	Start data download via PC-KIT-STICK
6		<i>[Cancel]</i>	Cancel data download via PC-KIT-STICK
7		<i>[Temperature warning min.]</i>	Set point lower temperature limit.
8		<i>[Target temperature]</i>	Set point operation temperature.
9		<i>[Min.]</i>	Showing minimum value of the temperature memory
10		<i>[Defrosting]</i>	Activate additional defrosting.
		<i>[Light]</i>	Switch interior lighting on/off permanently.
		<i>[Additional light / additional defrosting]</i>	Combination button: Activate additional defrosting. Switch interior lighting on/off permanently.



Depending on the version, certain buttons have a multiple function in combination with other buttons.

Tab. 2: Signs

Display	Description	Function
	'Defrosting'	Defrosting is active.
	'Humidity'	Humidity is active (temperature homogeneity improved, humidity high).
	'SuperFrost'	SuperFrost is active.

Display	Description	Function
	'Alarm'	An alarm was triggered.




2.2 Temperature sensor

The device is equipped with several temperature sensors.



PRO-ACTIVE control: The new sensors are PT1000 and not PTC100.

Please see the following table for the functions of the respective temperature sensors:

Colour	Designation	Function
	Evaporator sensor (red)	Controls the defrost phase
	Control sensor (brown)	Controls the monitoring unit
	Control sensor (violet)	Controls the monitoring unit

2.3 Error message and buzzer

Alarm

The unit is equipped with various alarm functions. When the alarm is triggered, it is optically and acoustically communicated via the display and control unit.

Temperature alarm or

- Occurs when the temperature exceeds or falls below the temperature warning limits.
- The remote warning contact is triggered.
- The display alternates between the current temperature and the error message.

Door Open alarm

- Happens if the door is open for more than 60 seconds.
- After 180 seconds, the door-open alarm is forwarded via the remote warning contact.
- The alarm goes off when the door is closed.

Power failure warning PF

- Happens when the unit is no longer being supplied with power.
- The display alternates between the current temperature and the error message.
- Units with a battery and PC-KIT record the temperature for approx. 30 hours after power is lost.
- Units with a battery record the min. and max. temperature for approx. 30 hours after power is lost.

Alarm indicating a defective display and control unit

- Happens when the display and control unit is no longer functioning.
- The unit must be taken out of operation.

Network-independent warning/ power loss



CAUTION!

Risk of damage to property and personal injury due to failure of the power supply.

The battery supplies power to the warning device only in the event of a power failure. The refrigeration unit is switched off, which may damage or destroy the chilled goods.

- After a power failure, ensure that the power supply to the refrigeration unit is restored.
- If necessary, move the chilled goods.

The unit has a battery that provides power to the warning unit for approx. 30 hours in the event of a loss of power. This allows the temperature memory to continue to record in the event of a loss of power and the temperature alarm is triggered outside of the temperature warning limits. Complete temperature documentation is guaranteed for every refrigerator and freezer.

The battery charges during normal operation and is monitored by an automatic charging system.






Replace the battery when the 'bATT' error message occurs. ↪ Chapter 4.11 'Replacing the battery' on page 37

2.4 Safety regulation

As long as the temperature is outside the set warning limits, the monitoring unit intervenes in temperature control. When the temperature is too low, the compressor circuit is interrupted using relay K1B.

2.5 Key switch

Tab. 3: Position of the key switch

Key position	Position	Function	Description
	"0"	Switch off unit.	<ul style="list-style-type: none"> ■ Switch off unit (for example to perform a restart). ■ The Display shows OFF. ■ To put the unit in a currentless state, unplug the power plug.
	"1"	Switch on unit.	<ul style="list-style-type: none"> ■ Operate unit. ■ Display and reset unit values (for example actual temperature or temperature limits). ■ Confirm alarms.
	"P"	Program unit.	<ul style="list-style-type: none"> ■ Adjust unit (for example target values or temperature warning limits). ■ Display parameter lists. ■ Confirm alarms.



Always operate unit in key position "1" in order to prevent manipulation. During normal operation, remove key and store it safely.

2.6 Defrosting

Automatic defrosting



The device defrosts automatically every 12 hours.

Automatic defrosting takes place with temporal and thermal monitoring.



When the defrosting process is active, the LED indicator 'Defrosting' lights up in the Display.

The defrost water is channelled into the external defrost water tank and evaporates there.

Additional defrosting

In addition to automatic defrosting, the defrosting process can be started manually.



- ➔ Hold down the [Additional defrosting] button for 4 seconds.
- ⇒ The defrosting process starts.

2.7 Data export on the device with PC-KIT-STICK

Saving data to USB

The refrigerator has the option of exporting the stored data to a standard USB stick (max. 32 GB, formatted in FAT32 format). A .jsn file is stored on the stick which can then be imported into KIRSCH-Datanet from version 5.0 by means of the procedure described below. The following procedure must be applied in order to save the data to the USB stick:










NOTICE!





The internal temperature memory records data for up to 90 days. It is therefore necessary to export the data regularly in order to avoid data gaps. We recommend exporting the data at regular intervals, e.g. weekly or monthly.



NOTICE!

Data can also be transmitted when the device is switched off (standby).

No:	Instruction	Display
1	Insert a USB stick	Wait until the following display appears:
		
		The controller is ready for data transfer (continue with 2.)
		The USB stick is full. Please clear it, or use a different stick.
		
		USB stick error when accessing the file system, or error when initialising the file system. Please use a different stick.
		
		USB stick error
		Data export failed, reinitialisation necessary. To do this, switch off the device using the key switch, then disconnect the device from the mains for 1 minute. Then reconnection the device and start data export
		
2	 Press the start button	The data is copied to the USB stick. ATTENTION: Never remove the USB stick during the copying process (see 5.). This will render the data otherwise unreadable/corrupt. If necessary, cancel the process first (see 4.).
		 

No:	Instruction	Display
3		The message "rdy" appears as soon as all data has been copied successfully. The USB stick can then be removed. 
4	 Cancelling with the Stop button	If desired, the copying process can be cancelled by pressing the Stop button. The USB stick can be removed as soon as "rdy" appears on the display. The data is then incomplete. 
5	Error message	If the stick is removed during copying, this message appears for approx. 5 seconds and the normal temperature display then reappears. 

Importing data into Kirsch Datanet

Establishing a connection

User name:

Password:

Server:

1. Open the KIRSCH Datanet software client and log in as Admin.

The default access details are:

- Username: Admin
- Password: admin



NOTICE!

Kirsch Service will not be able to reset the password if the default access details are changed and the password can no longer be found. You will then no longer have access to the database.

Refrigerator properties

Header data

Name:

Location:

☐ Monitoring on.

Acquisition rate: s

Connection

☒ USB-Gateway:

☐ IP address:

Refrigerator type

Refrigerator type:

Upper warning limit: °C Y-scale max: °C

Lower warning limit: °C Y-scale min: °C

Comment


2. Add a new refrigerator via the menu item "Verwaltung/Kühlschrank hinzufügen" [Administration/Add refrigerator].
3. The following window appears:
4. Enter the name of the device and, optionally, the location.
5. Uncheck the box next to "Monitoring on" [Monitoring on].
6. Set the USB export address if several devices are to be managed.



NOTICE!


Each address can be used only once.

7. Select the cabinet type.
8. Confirm by pressing "OK".




9. The newly created device appears in the "Kühlschränke" [Refrigerators] window on the left.

10. Insert the USB stick into the USB port on the computer.




11. Select the menu item "Datei/Datei importieren" [File/Import file] in the software.



12. Select the .json file to be imported from the USB stick via "Ordnersymbol" [Folder icon].

13. Select the desired refrigerator.

14. Click "Importieren" [Import] to import data.




15. The imported data can be displayed by double-clicking on the desired refrigerator.

3 Function check

3.1 Checking the alarm trigger

The temperature display is electronically delayed to correspond to the chilled goods temperature. The electronic delay must be switched off to test the alarm trigger.



1. ➔ Set the key switch to 'P'.



2. ➔ Press and hold [Temperature warning max.] and [Temperature warning min.] simultaneously for 4 seconds.

⇒ The display shows a flashing dot [.].

The test function has started.




For ProActive control, bear in mind that there is an upper (brown) and a lower (violet) sensor. Both sensors must be tested because the upper sensor warms up more quickly than the lower sensor. To test the upper sensor, press the 'Temperature warning max' button and to test the lower sensor, press the 'Temperature warning min' button.

3. ➔ Warm or cool the monitoring sensor in the interior to the appropriate temperature limit.

⇒ The buzzer sounds.

The remote warning contact is triggered.

The display flashes between the temperature and the error message.



4. ➔ Set the key switch to "1" to end the test function.

The test function ends automatically after 10 minutes.

⇒ The current temperature is shown on the display.

3.2 Checking sensors

1. ➔ Replace a defective sensor (e.g. F1) with any other sensor (e.g. F2).

- ⇒ ■ The error number of the replaced sensor (F2) is displayed: sensor (F1) is defective.
- No error is displayed: The sensor input on the control and monitoring unit is defective.

2. ➤ Replace the defective component.


➤ Chapter 4.5 'Replacing sensors' on page 35

3.3 Initialising the control and monitoring unit

By resetting the control and monitoring unit, Con, PAr and iEC errors can be resolved.



1. ➤ Set key switch to "0" to switch off the unit.
2. ➤ Unplug the power plug.
3. ➤ After 5 seconds, reconnect the power plug with the power grid.



4. ➤ Set the key switch to "1" to put the unit back into operation.
⇒ The internal memory of the control and monitoring unit is restarted.



If Con, PAr and iEC errors are triggered again within a short period of time, replace the control and monitoring unit (➤ Chapter 4.6 'Replacing the control and monitoring unit' on page 36).

3.4 Comparison of actual values

Carry out a comparison of actual values for the control in the following circumstances:

- When the measured interior temperature does not match the current display.
- When the measured interior temperature does not match the target value setting.
- When the control board or the control sensor and the display sensor have been replaced.


The interior temperature must be measured in the centre of the interior. When doing so, ensure that the sensor of the temperature measurement device is buffered so that the air temperature is not measured.

The control sensor and the display can be compared. The following rules apply for this:

- If the interior temperature is too warm, the control sensor setting must be represented as a negative number.
 - Example: Target temperature +4 °C, measured temperature +6 °C, comparison -2.0
- If the interior temperature is too cold, the control sensor setting must be represented as a positive number.
 - Example: Target temperature +4 °C, measured temperature +3 °C, comparison +1.0
- If the display shows a temperature that is too high, the display setting must be represented as a negative number.
 - Example: Displayed temperature +6 °C, measured temperature +4 °C, comparison -2.0
- If the display shows a temperature that is too low, the display setting must be represented as a positive number.
 - Example: Displayed temperature +3 °C, measured temperature +4 °C, comparison +1.0

If the measured temperature is different, first adjust the control sensor and then the display.

Selecting the Usr level parameter list



1. ➔ Set the key switch to 'P'.

2. ➔ Press and hold [Max.] and [Min.] simultaneously for 4 seconds.

⇒ The display shows **Rdr.**

3. ➔ Use [Max.] or [Min.] to select the Usr list.



4. ➔ Press [Target temperature] briefly.

⇒ The display shows **88.5**.



5. ➔ Use [Max.] to select one of the following parameters:

- **88.5**: Control sensor comparison
- **88.5**: Display comparison

Setting parameters




6. ➔ Press and hold [Target temperature].



7. ➔ Use [Max.] and [Min.] to change the value.



8. ➔ The value is saved by releasing [Target temperature].



9. ➤ Set the key switch to "1" to switch to operating mode.

Viewing correction of the offset value set



It is possible to view the offset correction set when the control sensor has been adjusted (parameter C66).

This is transferred to parameters H11 and H36 – therefore, parameter C66 is set to a value of 0.0 when it is viewed again.

This allows later verification of whether an adjustment has already been performed, and to what extent.

Proceed as follows to view the offset correction set:



1. ➤ Set the key switch to 'P'.



2. ➤ Press and hold [Max.] and [Min.] simultaneously for 4 seconds.

⇒ The display shows **Adr.**



3. ➤ Use [Max.] to navigate to point **PAE**.



4. ➤ Confirm [Target temperature].



5. ➤ Use [Max.] to navigate to level **H**.



6. ➤ Press and hold [Target temperature] and set **-12** using [Min.].



7. ➤ Use [Max.] to navigate to the preferred parameter, e.g. **H11**.




8. ➤ Press and hold [Target temperature] to view the offset set.

⇒ **In the background, the adjustment set under C66 is transferred automatically to the two relevant sensors H11 and H36.**



9. ➤ Release [Target temperature]. The set value is saved automatically.



10. ➤ Set the key switch to "1" to switch to operating mode.

4 Maintenance and Customization



DANGER!

Danger due to live parts!

- Unplug the power plug before carrying out any work.



WARNING!

Risk of death due to premature reactivation!

In the event of reactivation, there is a risk of death for anyone in the danger zone.

- Before reactivating, ensure that no one is in the danger zone.

Only carry out repair work with the power supply disconnected.

Protective equipment:

- Chemical-resistant protective gloves
- Protective gloves
- Protective goggles

4.1 Status displays and error messages




4.1.1 Status displays


Status displays are indications designed to provide information to the user (e.g. of ongoing defrosting).

No audible signal (hereinafter referred to as "buzzer") is emitted when the status is displayed.

A status display does not require immediate action by the user.

Tab. 4: Status displays

Display	Buzzer	Description	Measure	Unit key required
	–	LED display <i>[Defrosting]</i> lights up: ■ The defrosting process is active.	–	–
	✓	LED display <i>[Alarm]</i> lights up: ■ One or more alarms have been triggered (collective alarm).	–	–
	–	LED display <i>[SuperFrost]</i> lights up: ■ Super Frost is active.	–	–

Display	Buzzer	Description	Measure	Unit key required
	–	Standby display: <ul style="list-style-type: none"> The unit is connected to the power grid and the key switch is set to position "0". 	Turn the key switch to position "1" to switch on the device.	✓

4.1.2 Error messages



Error messages indicate a device malfunction.

Error messages alternate on the display with the temperature display.

If there is more than one error, the errors are shown one after the other on the display.





In addition to the indication on the display, the buzzer sounds to indicate the error.

The device indicates the following errors visually and audibly:

- Unit error ( 'Error messages from the device' on page 29)
- Software errors ( 'Error messages from the device' on page 29)

If error messages occur, proceed as described below:

Procedure for error messages

-  Deactivate the buzzer.
-  Assess the error display according to the table.
-  Carry out recommended measures.
-  Acknowledge the alarm message.



NOTICE!

Danger to chilled goods due to defective or faulty unit!

A defect or fault in the unit means that its cooling performance is no longer ensured. Reduced cooling performance can cause considerable damage to chilled goods.

- Select an alternative storage location for the chilled goods.
- Ensure operating and storage conditions.
- Transfer chilled goods to new location.



CAUTION!

Danger of death due to interrupted cooling chain during stock transfer!

The unit is used for commercial storage of pharmaceuticals, vaccines and pharmaceutical ingredients that require cooling in line with the specifications of the manufacturer. If the cooling chain of the chilled goods is interrupted by a stock transfer, the prescribed storage conditions are no longer met.

- Do not subject chilled goods to light during the stock transfer.
- Do not place chilled goods near radiators during the stock transfer.
- Make sure that chilled goods are stored in the replacement unit according to the specifications of the relevant manufacturer.



NOTICE!

Meaning of "X" for error and status messages



X is not shown on the display.


- Instead, the display shows a number that describes the relevant part.














In normal operation, USB sticks may not be connected to the unit because error displays and acoustic warning signals then do not occur.

Tab. 5: Error messages from the device






Display	Buzzer	Description	Measure	Unit key required
	✓	Sensor X: ■ Error or short circuit of the corresponding sensor. ■ The cooling controller is running the emergency program.	■ Check the sensors and connections at the control board. – Switch the sensors at the control board to test them. – If error message F1L is displayed, switch sensors F1 and F2 at the control board. – If error message F2L does not appear, replace the control board.	✓
	✓	Sensor X: ■ Error or failure of the corresponding sensor. ■ Cooling controller is running the emergency program.		✓


Display	Buzzer	Description	Measure	Unit key required
F 1 1 H	✓	<ul style="list-style-type: none"> Defective battery connection cable. Battery monitoring is not guaranteed. 	<ul style="list-style-type: none"> Switch off the alarm. Replacing the battery connection cable Replace the battery ( Chapter 4.11 'Replacing the battery' on page 37). 	✓
L X L	✓	Fan X: <ul style="list-style-type: none"> Speed of corresponding fan too low. The temperature of the chilled goods may fluctuate. 	<ul style="list-style-type: none"> Check the interior blower for resistance from the outside. If there is no blockage, replace the interior blower. 	✓
L X H	✓	Fan X: <ul style="list-style-type: none"> Speed of corresponding fan too high. The temperature of the chilled goods may fluctuate. 		✓
F A 1	✓	Fan: <ul style="list-style-type: none"> Fan does not reach the required minimum speed after restarting the device. The temperature of the chilled goods may fluctuate. 		✓
d F A	✓	Fan: <ul style="list-style-type: none"> Fan speeds deviate too much from each other. The temperature of the chilled goods may fluctuate. 		✓
r 0 1 L	✓	<ul style="list-style-type: none"> Measured power consumption of cooling machine too low. The temperature of the chilled goods may fluctuate. 	<ul style="list-style-type: none"> Carry out a control reset. Check the functioning of the control board and cooling machine if the error occurs again after a short time. 	✓
r 0 1 H	✓	<ul style="list-style-type: none"> Measured power consumption of cooling machine too high. The temperature of the chilled goods may fluctuate. 		✓
P A r	✓	Synchronisation error: <ul style="list-style-type: none"> Synchronisation error between control panel and monitoring circuit. No reliable cooling controller function. 	<ul style="list-style-type: none"> Carry out a control reset. Replace the control board and display if the error occurs again after a short time. 	✓

Display	Buzzer	Description	Measure	Unit key required
	✓	Connection problem: <ul style="list-style-type: none"> ■ Synchronisation error between control panel and monitoring circuit. ■ No reliable cooling controller function. 		✓
	✓	Connection problem: <ul style="list-style-type: none"> ■ Synchronisation error between control and alarm circuit. ■ No reliable cooling controller function. 		
	✓	Control error: <ul style="list-style-type: none"> ■ Error during self-testing in cooling controller. ■ The monitoring circuit is responsible for temperature control. 		✓
	✓	Door Open alarm: <ul style="list-style-type: none"> ■ The door has been open for more than 60 seconds. ■ The door contact switch is defective. 	<ul style="list-style-type: none"> ■ Close the door. ■ Check the function of the door contact switch. ■ Replace the door contact switch ( Chapter 4.10 'Replacing the door contact' on page 37). 	–
	✓	Defective battery: <ul style="list-style-type: none"> ■ Battery must be replaced. ■ The temperature documentation and alarm fail in the event of power failure. 	<ul style="list-style-type: none"> ■ Switch off the alarm. ■ Replace the battery ( Chapter 4.11 'Replacing the battery' on page 37). 	–
	✓	Power failure: <ul style="list-style-type: none"> ■ Device does not cool. ■ Alarm is active. ■ The remote warning contact is triggered. 	<ul style="list-style-type: none"> ■ Move chilled goods. ■ Check the power supply. 	–
	✓	Temperature alarm (high): <ul style="list-style-type: none"> ■ The upper temperature warning limit is reached or exceeded (for example due to very warm chilled goods or because the door was open too long). 	<ul style="list-style-type: none"> ■ Move chilled goods. ■ View and check the temperature warning limit. ■ If necessary, correct the values of the temperature warning limits. ■ Observe the temperature curve. 	✓

Display	Buzzer	Description	Measure	Unit key required
	✓	<p>Temperature alarm (low)</p> <ul style="list-style-type: none"> ■ The lower temperature warning limit has been reached or exceeded (e.g. after the door has been open for a long time when the cooling machine is cooling the device). ■ Safety device responds, monitoring circuit has switched off cooling machine. 	<ul style="list-style-type: none"> ■ Move chilled goods. ■ View and check the temperature warning limit. ■ If necessary, correct the values of the temperature warning limits. ■ Observe the temperature curve. 	✓
	✓	<p>Incorrect time during data output:</p> <ul style="list-style-type: none"> ■ In the event of an rtc error, the date and the year must be set. 	<ul style="list-style-type: none"> ■ Set the time, date and year. 	✓

Tab. 6: Error messages from the device

Display	Buzzer	Description	Measure	Unit key required
	–	<p>Status display:</p> <ul style="list-style-type: none"> ■ Copying process is running. ■ XX stands for the percentage progress of the copying process. 	<ul style="list-style-type: none"> ■ Do not remove the USB stick. 	–
	–	<p>Status display:</p> <ul style="list-style-type: none"> ■ The copying process is complete. 	<ul style="list-style-type: none"> ■ Remove the USB stick. 	–
	–	<p>Status display:</p> <ul style="list-style-type: none"> ■ The memory of the USB stick is full. 	<ul style="list-style-type: none"> ■ Use a USB stick with sufficient storage capacity. 	–
	–	<p>Error message:</p> <ul style="list-style-type: none"> ■ USB stick error when accessing the file system, or error when initialising the file system. 	<ul style="list-style-type: none"> ■ Please use a different stick 	–
	–	<p>Error message:</p> <ul style="list-style-type: none"> ■ Data export failed, reinitialisation necessary. 	<ul style="list-style-type: none"> ■ Switch off the device using the key switch, then disconnect the device from the mains for 1 minute. ■ Then recommission the device and start data export. 	–


Display	Buzzer	Description	Measure	Unit key required
	—	Error message: ■ USB stick was removed during the copying process.	■ Remove the USB stick and try again after 1 minute.	—



NOTICE!

Note that the USB stick must be removed after data output as otherwise functions may be restricted.

4.1.3 Setting the time and date



1. ➔ Set the key switch to 'P'.



2. ➔ Press the [Target temperature] button twice.

⇒ In 24-hour mode, the display shows the time for approx. 3 seconds.



3. ➔ Use the [Min.] or [Max.] button to set the desired time.

4. ➔ Confirm using the [Target temperature] button.



Press the [Target temperature] button several times until the desired setting, such as e.g. time, year or date, appears on the display. After the desired change has been made using the [Min.] or [Max.] button, confirm the entry using the [Target temperature] button.

The following setting options are available:

- Time (24-hour mode)
- Year (format: YYYY)
- Date (format: DD.MM)

⇒ The time, year or date is set.

5. ➔ Set the key switch to "1".

4.2 Switching the door hinge

Changing the door hinge from right to left

The door must be turned 180° to ensure that the installed spring compressors still function correctly after the hinge switch.

1. ➤ Unscrew the machine compartment panel and move it back.
2. ➤ Open the door to 90°.
3. ➤ Undo the upper hinge.
4. ➤ Tip the door forward until the hinge releases out of the opening.
5. ➤ Lift the door out of the lower hinge.
6. ➤ Unscrew the lower hinge.
7. ➤ Screw the previously upper hinge into the position for the lower hinge on the other side.
8. ➤ Turn the door 180° and with an opening angle of 90°, place it onto the hinge that is now in the lower position.
9. ➤ Guide the other hinge into the housing opening.
10. ➤ Mount the upper hinge.
11. ➤ Close the door and ensure it is aligned flush.
12. ➤ Switch the door handle as described in [Chapter 4.3](#) 'Switching the door handle' on page 35.

Changing the door hinge from left to right

1. ➤ Unscrew the machine compartment panel and move it back.
2. ➤ Undo the upper hinge.
3. ➤ Tip the door forward until the hinge releases out of the opening.
4. ➤ Lift the door out of the lower hinge.
5. ➤ Unscrew the lower hinge.
6. ➤ Screw the previously upper hinge into the position for the lower hinge on the other side.
7. ➤ Turn the door 180° and place it onto the hinge that is now in the lower position.
8. ➤ Manually pretension the other hinge and guide it into the housing opening.
9. ➤ Mount the upper hinge.
10. ➤ Close the door and ensure it is aligned flush.
11. ➤ Switch the door handle as described in [Chapter 4.3](#) 'Switching the door handle' on page 35.

4.3 Switching the door handle

1. ➤ Remove the two white plugs and the insulator on the inside of the door.
2. ➤ Unscrew the fixing screws.
3. ➤ Unscrew the threaded bushing.
4. ➤ Remove the door handle.
5. ➤ Insert the new door handle into the door from the outside.
6. ➤ Screw in the threaded bushing but do not tighten it. When doing so, ensure that no foam material enters the plunger guide.
7. ➤ Screw in the fixing screws.
8. ➤ Insert the isolator and plugs.
9. ➤ Tighten the threaded bushing slightly.

4.4 Replacing the display unit

Display unit above the door in the machine compartment panel

Depending on the model, the display unit is in the panel above the door or lower down in the machine compartment panel.

1. ➤ Unplug the connection plug for the control and monitoring unit.
2. ➤ Undo the screws of the display holder.
3. ➤ Remove the display unit from the machine compartment panel.
4. ➤ Fit the new display unit in the reverse order.

4.5 Replacing sensors

Guide the defective sensor out of the housing

Access to the control and monitoring unit behind the cover plate behind the cooling machine

Detach the defective sensor from the control and monitoring unit

Mounting a new sensor

1. ➤ Undo the fixing screws of the defective sensor.
2. ➤ Remove the sealing compound at the opening for the sensor.
3. ➤ Guide the defective sensor out of the housing.
4. ➤ Undo the cover plate screws behind the cooling machine.
5. ➤ Remove the cover plate.
6. ➤ Remove the connection plug of the defective sensor from the board.
7. ➤ Mount the new sensor in the reverse order. When doing so, ensure that the opening for the sensor is closed tightly with a new sealing compound.

4.6 Replacing the control and monitoring unit

Depending on the model, the control and monitoring unit is under the cover plate next to/behind the cooling machine or in the lower part of the unit behind the fitting cover.

Control and monitoring unit behind the cover plate at the side of the cooling machine

1. ➤ Undo the cover plate screws next to the cooling machine.
2. ➤ Remove the cover plate.
3. ➤ Remove the connection plug to the display unit.
4. ➤ Undo the screws of the control and monitoring unit on the terminal board.
5. ➤ Remove the control and monitoring unit.
6. ➤ Fit the new control and monitoring unit in the reverse order.
7. ➤ Recalibrate parameters 'C66' and 'J11' ➤ *Chapter 3.4 'Comparison of actual values' on page 24.*

4.7 Replacing the interior blower

1. ➤ Undo the fixing screws of the ventilation plate.
2. ➤ Remove the ventilation plate.
3. ➤ Remove the plug connectors between the interior blower and the connection cable.
4. ➤ Detach the interior blower and the mounting plate from the back wall and remove.
5. ➤ Remove the interior blower from the mounting plate.
6. ➤ Secure the new interior blower to the mounting plate.
7. ➤ Reinstall the new interior blower on the mounting plate.

4.8 Replacing the door lock

Door lock above the door in the top section

1. ➤ Open the door.
2. ➤ Undo and remove the 2 screws beside the deadbolt.
3. ➤ Undo the screw in the machine compartment panel and remove it.
4. ➤ Remove the lock.
5. ➤ Remove the mounting bracket.
6. ➤ Fit the new lock in the reverse order.

4.9 Replacing the door seal



The door does not have to be unscrewed to replace the door seal.

1. ➤ Remove the defective door seal from the groove.
2. ➤ Clean the groove of any residue.
3. ➤ Depending on the model, spray fresh silicon into the groove.
4. ➤ Insert the new seal into the groove and gently tap it into the groove (e.g. with the ball of your thumb or with a mallet).
5. ➤ Ensure that the seal sits in the groove to the cover lip and seals completely at the housing.

4.10 Replacing the door contact

Magnetic door contact switch in the machine compartment panel above the door (reed contact)

1. ➤ Unplug the plug connector.
2. ➤ Open the door.
3. ➤ Push the clamp springs back with the screwdriver.
4. ➤ Pull out the switch at the pin.
5. ➤ Insert a new switch.
6. ➤ Secure the new switch with clamp springs.

4.11 Replacing the battery

Battery at the top of the top unit next to the cooling machine

1. ➤ Disconnect the red cable +.
2. ➤ Remove the red and blue cable from the battery.
3. ➤ Pull the battery out of the holder.
4. ➤ Install the replacement battery in reverse order.

4.12 Replacing the cold room lighting (optional)

The cold room lighting consists of a LED light rail. It can only be replaced as a complete unit including the lighting ballast. The lighting ballast is attached at the control unit.

1. ➤ Undo the cover plate screws behind the cooling machine.
2. ➤ Remove the cover plate.
3. ➤ Release the metal bracket at the ends of the light strip from the unit.

4. ➤ Remove the metal bracket.
5. ➤ Remove the lighting ballast.
6. ➤ Remove the sealing compound on the cable bushing.
7. ➤ Guide the cable out of the housing.
8. ➤ Fit the new LED light strip in the reverse order.

4.13 Replacing the cooling machine



WARNING!

Danger due to unauthorised intervention in the cooling circuit!

The cooling machine contains the natural refrigerant propane R290 / isobutane R600a. The refrigerant used in the unit is flammable and can cause serious injury if it comes into contact with the eyes or skin. Unauthorised intervention in the cooling circuit poses a danger of injury!

- On site, only the entire cooling machine can be replaced.

The cooling machine is a permanently technically sealed system in line with EN 1127-1.

For Eco units: Only the entire cooling machine can be replaced. The refrigerant R600a or R290 is used in the cooling circuit of the Eco unit.

1. ➤ Remove the machine compartment panel.
2. ➤ Open the door.
3. ➤ Remove the ventilation plate.
4. ➤ Remove the upper right-hand line cover in the interior.
5. ➤ Remove the upper cover section.
6. ➤ Open the upper right-hand line opening.
7. ➤ Unscrew the compressor set.
8. ➤ Remove the machine cable at the terminal box of the cooling machine.
9. ➤ Release the evaporator and remove it to the front.
10. ➤ Carefully remove the cooling machine to the front. When doing so, ensure that the tubing is not damaged or kinked.



For easier installation, the cold lines can be fixed to the compressor set at the condenser.

- 11.** ➤ Fit the new cooling machine in the reverse order. When doing so, ensure that the tubing is not kinked or damaged.

Replacing the cooling machine

5 Customizing parameters




NOTICE!

Material damage from incorrectly set parameters!

The basic settings of the unit can be changed using the following parameter lists. This could result in severe damage to the chilled goods or the unit.

The parameter settings may only be carried out by trained staff.

Calling up the parameter list



1. ➔ Set the key switch to 'P'.



2. ➔ Press and hold [Max.] and [Min.] simultaneously for 4 seconds.

⇒ The display shows **Rdr**.



3. ➔ Use [Max.] or [Min.] to navigate to "Usr" level.



4. ➔ Press [Target temperature] to confirm the level.



5. ➔ For entries in "Usr", [Max.] and [Min.] can be used to navigate through the parameter list.

Setting parameters



6. ➔ Select the parameter to be changed.

7. ➔ Press and hold [Target temperature].

⇒ The currently set value for the parameter is set.



8. ➔ Set the parameter value using [Max.] and [Min.].




9. ➔ Release [Target temperature].

⇒ The set value is saved.

The display changes automatically back to the basic state if no buttons are pressed for 45 seconds.

Set the key switch to "1" after finishing work.



Usr level

Parameter	Function description	Adjustment range	Factory setting
H65	Display of power consumption	Read only	—
Rn1	Cold room sensor (Mixed sensor delayed)	Read only	—
Rn2	Superfrost sensor (not used)	Read only	—
Rn3	Moisture sensor (not used)	Read only	—
Rn4	Cold room sensor, upper	Read only	—
Rn5	Cold room sensor, lower	Read only	—
Rn6	Evaporator sensor	Read only	—
Rn7	T-limit sensor alarm, upper	Read only	—
Rn8	T-limit sensor alarm, lower	Read only	—
Rn9	T-prevention sensor, upper	Read only	—
RnA	T-prevention sensor, lower	Read only	—
C51	Light function	0: off 1: on	0
C52	Light delay time	0...300 s	10 s
C64	Compressor shut-down hysteresis	-5.0 K...5.0 K	0.0
C66	Control sensor comparison	-5.0 K...5.0 K	0.0
F8	Fan speed regular operation	0...10000	2100
b5	Battery buffering	0: Deactivated 1: Always active 2: Active when system on 3: Active for key switch 1 or 2	3

Parameter	Function description	Adjustment range	Factory setting
d1	Defrosting interval	0: no automatic defrosting 1...99 h	12
d3	Stop at defrost temperature	0.0...30.0 °C	5.0
d7	Temperature difference for previous cooling	-15,0...0,0 K	-0.5
d11	Drip tray heating delay time	0...60 min	10
R25	"PowerFail" alarm suppression time	0...300 s	15
G10	Recording interval for data logger	0...900 s	120
G95	USB export ID of the refrigerator	0...9999	0
r5	Automatic summer/winter time	0: Deactivated 1: Europe 2: USA 3: Australia 4: New Zealand	1
r6	Time zone (UTC) [x * 15 min]	-48...52 (-12 h...13 h)	4
J11	Display comparison	-20.0...20.0 K	0.0
L48	TCP/IP activation	0: deactivated 1: activated	0
L50	IP address 3	0...255	192
L51	IP address 2	0...255	168
L52	IP address 1	0...255	0
L53	IP address 0	0...255	101
L55	IP mask 3	0...255	255
L56	IP mask 2	0...255	255
L57	IP mask 1	0...255	255

Parameter	Function description	Adjustment range	Factory setting
L58	IP mask 0	0...255	0
L60	IP gateway 3	0...255	192
L61	IP gateway 2	0...255	168
L62	IP gateway 1	0...255	0
L63	IP gateway 0	0...255	200
L70	MAC address 5	0...0xFF Read only	—
L71	MAC address 4	0...0xFF Read only	—
L72	MAC address 3	0...0xFF Read only	—
L73	MAC address 2	0...0xFF Read only	—
L74	MAC address 1	0...0xFF Read only	—
L75	MAC address 0	0...0xFF Read only	—

6 Appendix

The appendix contains the following applicable documents:

- Technical data
 - ↳ Chapter 6.1 'Technical data' on page 46
- wiring diagram
 - ↳ Chapter 6.2 'wiring diagram' on page 49
- Exploded view
 - ↳ Chapter 6.3 'Exploded view' on page 51
- Spare parts list
 - ↳ Chapter 6.4 'Spare parts list' on page 52

6.1 Technical data

	BL 100 PRO-ACTIVE	BL 176 PRO-ACTIVE	BL 300 PRO-ACTIVE/ ULTIMATE*	BL 520 PRO-ACTIVE/ ULTIMATE*
Cooling capacity, litres	90	170	280	500
Temperature setting (approx.), °C	+4	+4	+4	+4
Voltage, V	220 – 240	220 – 240	220 – 240	220 – 240
Frequency, Hz	50	50	50	50
Refrigerant quantity, grams	42	30	30	30
Refrigerant	R600a	R600a	R600a	R600a
Current consumption, amps	0.42	0.53	0.43	1.53
Power consumption, watts	92	120	95	250
Normal consumption, kWh/24	0,64	1.22	0.63	1.50
Perm. ambient temperature, °C	+10 to +38	+10 to +38	+10 to +38	+10 to +38
External dimensions incl. distance to wall (WxDxH), cm	54 x 54 x 82	70 x 77 x 122	68 x 72 x 131.5	77 x 76 x 195.5
Usable sizes (WxDxH), cm	39.5 x 34 x 34	51 x 34.5 x 46	49 x 39 x 92 (Usable depth 15 cm less at the bottom)	56 x 39.5 x 10
External dimensions with door open 90° (WxD), cm	54 x 106	70 x 136	67 x 133	77 x 144
Shelf size (WxD), cm	–	–	–	–
Internal drawer dimensions (WxDxH), cm	39.5 x 34 x 10	51 x 34.5 x 10	49 x 39 x 10 top drawers 51 x 31 x 10 bottom drawer	56 x 39 x 10
Max. drawer/shelf load, kg	25 / -	50 / -	50 / -	50 / -
Weight net/gross, kg	48 / 57	77 / 87	81 / 93	138 / 161
Noise emission, dB(A)	40.5	38.6	40.3	41.2

	BL 720 PRO-ACTIVE/ ULTIMATE*	FROSTER- BL-178	FROSTER- BL-180	FR BL 330 PRO-ACTIVE/ ULTIMATE*
Cooling capacity, litres	700	170	170	300
Temperature setting (approx.), °C	+4	-32	-42	approx. -32
Voltage, V	220 – 240	220 – 240	220 – 240	220 – 240
Frequency, Hz	50	50	50	50
Refrigerant quantity, grams	90	400	480	440
Refrigerant	R600a	R404a	R404a	R404a
Current consumption, amps	1.53	3.00	2.26	3.2
Power consumption, watts	250	500	453	550
Normal consumption, kWh/24	1.56	4.88	6.54	6.59
Perm. ambient temperature, °C	+10 to +38	+10 to +32	+10 to +32	+10 to +32
External dimensions incl. distance to wall (WxDxH), cm	77 x 98 x 195.5	70 x 77 x 122	70 x 77 x 122	74 x 77 x 159 74 x 77 x 166*
Usable sizes (WxDxH), cm	56 x 60.5 x 129	55 x 52 x 60 Usable depth 11 cm, Usable width 3 cm, Usable height 13 cm	55 x 52 x 60 Usable depth 15 cm, Usable width 3 cm, Usable height 13 cm	54 x 37 x 80
External dimensions with door open 90° (WxD), cm	77 x 166	70 x 136	70 x 136	74 x 142
Shelf size (WxD), cm	–	–	–	–
Internal drawer dimensions (WxDxH), cm	56 x 60.5 x 10	50 x 34 x 10	50 x 32 x 10	54 x 37 x 10
Max. drawer/shelf load, kg	50 / -	50 / -	50 / -	50 / -
Weight net/gross, kg	164 / 190	92 / 105	100 / 110	120 / 135
Noise emission, dB(A)	41.8	50.3	51.2	50.4

	FR BL 530 PRO-ACTIVE/	FROSTER- BL-650	FR BL 730 PRO-ACTIVE/
Cooling capacity, litres	500	650	700
Temperature setting (approx.), °C	approx. -32	-41	approx. -32
Voltage, V	220 – 240	220 – 240	220 – 240
Frequency, Hz	50	50	50/60
Refrigerant quantity, grams	620	820	620
Refrigerant	R404a	R404a	R404a
Current consumption, amps	3.8	6.5	4
Power consumption, watts	740	990	790
Normal consumption, kWh/24	10.81	11.10	11.1
Perm. ambient temperature, °C	+10 to +32	+10 to +32	+10 to +32
External dimensions incl. distance to wall (WxDxH), cm	77 x 76 x 195.5	83 x 102 x 197	77 x 97 x 195.5
Usable sizes (WxDxH), cm	56 x 39.5 x 129	62 x 76 x 131 Usable depth 26 cm, Usable width 21 cm, Usable height 14 cm lower	56 x 60.5 x 129
External dimensions with door open 90° (WxD), cm	77 x 144	86 x 177	77 x 166
Shelf size (WxD), cm	–	–	–
Internal drawer dimensions (WxDxH), cm	56 x 39.5 x 10	59 x 56 x 22 top 59 x 62 x 22 bottom	56 x 60.5 x 10
Max. drawer/shelf load, kg	50 / -	50 / -	50 / -
Weight net/gross, kg	168 / 191	220 / 260	188 / 214
Noise emission, dB(A)	51.5	62.3	49.7

* Details to follow

** Details not available (other standard)

6.2 wiring diagram




Fig. 3: 500-022_Wiring_Diagram_BL_520_720_PRO-ACTIVE_CR



The wiring diagram configuration varies depending on the model. The legend describes the maximum configuration.

A1	Operating and display unit	T4	Power supply (230V AA/ 10VDC)
A2	Control and monitoring board	V1	Drip tray heating
B1	Safety barrier 9001 01-086-150-10	V2	Melt water evaporation
B2	Safety barrier 9001 01-086-150-10	V3	Interior lighting (option)
B3	Safety barrier 9001 01-086-150-10	V4	Door frame heating
B4	Pressure switch	W1	Remote warning contact connection
B5	Pressure switch (pressure side)	W2	Battery connection
B6	Pressure switch (suction side)	W3/W4	Connection for power, compressor, ...
B7	Temperature limiter	W5	Drip tray heating connection
C1	Condenser	W6/W7	DC fan connection
C2	Capacitor for cooling machine fan/motor operating capacitor	W8/W9	Connection for sensor, door contact
F1	Temperature sensor control / monitoring	W10	Internal bus
F2	Evaporator temperature sensor	W11	LAN connection
F6	Temperature sensor control / monitoring	W12	Display connection
G1	Battery (option)	X1	Terminal for FU signal
K1A	Compressor relay (hybrid)	X3	Cooling machine terminal
K1B	Compressor relay (hybrid)	X10	Remote warning contact
K2	Defrosting relay	X13	Power plug
K3	Melt water evaporation relay, door frame heating	X15	Cooling machine terminal
K4	AC interior fan relay	X16	3-way socket
K5	Lighting relay	Y1	Magnetic valve for hot gas line
K6	Relay for drip tray heater/condenser	Y2	Magnetic valve
K7	Remote maintenance relay	Y3	Vacuum compensation valve
K10	Start-up relay	1F1	Residual-current circuit breaker
K11	Finder relay 30A		
M1	Compressor motor	brown (braun)	
M2	Condenser blower	blue (blau)	
M3	Recirculation blower	black (schwarz)	
M4	Recirculation blower	grey (grau)	
S1	Door contact	green-yellow (grün-gelb)	
T1	Power supply (230V AA/ 10VDC)	white (weiß)	
T2	Lighting ballast (option)	red (rot)	
T3	Power supply (230V AA/ 10VDC)	purple (lila)	

6.3 Exploded view




Fig. 4: 700-052-2_Exploded_View_BL 720-ECO PRO-ACTIVE also_CR

6.4 Spare parts list

Cooling machine

Current item number	Meaning
10283	Cooling machine VKD 4703 KS
10961	Relay, Klixon, terminal block
10960	Start capacitor
10964	Blower motor
10963	Metal protective grating
10965	Capacitor fan
114660	Rubber buffer

Control

Current item number	Meaning
500319	Control and monitoring unit (program J1 = 19)
500369	Display unit (program J1 = 19)
11441	Temperature sensor F1, 2 m, brown (not depicted)
11439	Temperature sensor F2, 2 m, red (not depicted)
11442	Temperature sensor F6, 3 m, violet (not depicted)
114391	Main body, 3-pole for remote warning contact
114392	3-pin plug for remote warning contact
114407	Key for key switch (not depicted)
202067	6 V 4 Ah battery
100413	Battery bracket
500489	Battery holder/interface angle
202084	Magnetic switch (reed contact)
12025	Evaporation heating

Housing

Current item number	Meaning
500435	Machine compartment panel
500465	Stainless steel machine compartment panel
254138	Fitting cover with filmed front
500436	Lock fastener (not depicted)
253142	Lock for door, closing 1001
270125	Key no. 1001
500487	Terminal board
500408	Terminal board cover plate
500407	Melt water evaporation dish
253123	Hinge with square
253122	Hinge with pin Ø = 7 mm
114113	Fixed castor
113164	Castor

Interior

Current item number	Meaning
11369	Axial blower 24 Volt DC
500486	Cross plate for blower
500430	Upper ventilation plate
500488	Lower ventilation plate
500448	Left attachment rails for ventilation plate
500449	Right attachment rails for ventilation plate
253343	Support rails 1200 mm
270106T1	Drip tray 480 x 85 x 50/30 mm, without heating cable
224624	Angle bracket for door handle
700402	Cover strip
11832	Drawer with differential pull-out mechanism
271628LI	Differential pull-out mechanism with support rail, left
271628RE	Differential pull-out mechanism with support rail, right (not depicted)
11851	Roller drawer guides, MS, left, differential pull-out mechanism 600 mm.
11852	Roller drawer guides, MS, right, differential pull-out mechanism 600 mm (not depicted)
500402	Line duct cover part

Refrigerator door

Current item number	Meaning
254027	Door, white, right/left, 1539 x 755 mm, self-closing
254028	Door, stainless steel, right/left, 1539 x 755 mm self-closing
253139	Clamp for door
270113	Door seal 1496 x 716 mm
114586	Door handle 50 mm
253124	Plastic door bearing bush
270141	Plastic bushing for hinge pin
202082	Magnet for door contact switch
202037	Energy storage unit (not depicted)

Option (not depicted)

Current item number	Meaning
10421	Glass door, right/left, 1539 x 755 mm, self-closing
10556	Manual emptying, 720 ltr.
10430	LED light strip, complete
11947	PC-KIT-NET
253217	Melt water dish for manual emptying
271193	Drawer divider set consisting of longitudinal divider 602 mm, cross divider 100 mm and locking strip 490 mm, model BL-720
271191	Longitudinal divider 602 mm, model BL-720
114624	Locking strip, 490 mm
114595	Cross divider with card pocket 100 mm long

7 Index

A

Additional defrosting	19
Appendix	45
Applicable documents	45

B

Button functions (overview)	15
Buzzer	17

C

Changing the door hinge	34
Character	8
Check alarm triggering	23
Check functions	
Actual value	24
Alarm triggering	23
Checking sensors	23
Comparison of actual values	24
Control and monitoring unit	
replace	36
Resetting	24

D

Data export	20
Defrosting (automatic)	19
Defrosting (manual)	19
Display	15
Display adjustment	24

E

Error messages	17
Error messages (overview)	28
Exploded view	51

F

Foreseeable misuse	8
Function check	
Control and monitoring unit	24
Sensors	23

I

Images on the display (overview)	15
Intended use	8

K

Key switch	15, 19
----------------------	--------

M

Maintenance and Customization	27
Door seal	37

O

Operation	15
---------------------	----

P

Parameter list	42
Parameterisation	41
Personnel	12
Protective equipment	13
Protective gloves	13

R

Rating plate	11
Repair	
Battery	37
Cold room lighting (optional)	37
Control and monitoring unit	36
Cooling machine	38
Display unit	35
Door contact	37
Door handle	35
Door lock	36
Interior blower	36
Sensors	35
Replacing cold room lighting (optional)	37
Replacing sensors	35
Replacing the battery	37
Replacing the cooling machine	38
Replacing the display unit	35
Replacing the door contact	37

Replacing the door handle (optional)	35
Replacing the door lock	36
Replacing the door seal	37
Replacing the interior blower	36
Residual risks	8

S

Safety	7
Safety regulation	18
Spare parts list	51
Control	52
Cooling machine	52
Housing	53
Interior	54
Optional	55
Refrigerator door	54
Status displays (overview)	27
Symbols	7

T

Technical data	46
Temperature sensor (overview)	17

U

USB	
Data import	21
USB port	15
Usr	42

W

Warning sound	17
Wiring diagram	49